

**STOPS BLEEDING  
IN SECONDS**



**VETIGEL**<sup>®</sup> A fast-working hemostatic gel for use in veterinary procedures

**Benefits for clinics**

- Reduces costly O.R. times
- Speeds up recovery
- Decreases the risks to your patients with pre-existing conditions
- Faster patient turnaround time
- Lower cost of hemostasis

**Benefits for patients**

- Reduces blood loss
- Shortens time under anesthesia
- Lower risk of systemic shock
- Avoids risk of post-op bleeding
- Being reunited faster with pet owner

**VETIGEL**<sup>®</sup> Products



**VETIGEL**<sup>®</sup>  
Hemostatic Gel  
(5 Pack)



**VETIGEL**<sup>®</sup>  
Hemostatic Gel Kit  
(2 Pack with  
Accessories)




**VETIGEL**<sup>®</sup>  
Accessories  
(4 Dental Tips,  
2 Adapters)


**Where can VETIGEL**<sup>®</sup> be used?

**VETIGEL**<sup>®</sup> can be used on a variety of procedures, including but not limited to: dental extractions, oral mass removals, abdominal and thoracic mass removals, soft organ biopsies, enucleations, spays & neuters, amputations, venous & arterial bleeds, cholecystectomies, catheterization, lacerations, and ear & tail docking.


**DENTAL EXTRACTIONS**  
(10-15 TEETH)



**AMPUTATIONS**



**ENUCLEATIONS**



\*Study Time Savings vs. Standard of Care  
Dental extractions: 15 min. | Amputations: 20 min. | Enucleations: 5 min.

# Frequently Asked Questions

## How does VETIGEL® work?

VETIGEL® is a hemostatic gel that immediately stops bleeding via mechanical action, without the need for applied pressure. Gel that comes directly in contact with blood ionically crosslinks, forming a strong mechanical barrier that maintains durable and long-term hemostasis at the wound site.

VETIGEL® allows the patient to rapidly produce their own stable endogenous fibrin patch at the wound site. VETIGEL® can be easily removed without disturbing the fibrin patch thanks to its flowable, non-porous nature.

## Is VETIGEL® absorbable?

VETIGEL® is absorbable, and is intended to be absorbed within weeks of application, while minor residuals may take longer to absorb.

## What is VETIGEL® made of?

VETIGEL® is a patented flowable hemostatic gel comprised of algae- and fungi-derived biopolymers. Its components are sodium alginate, poly (N-acetyl-D-glucosamine, D-glucosamine) and water.

## How quickly does VETIGEL® stop bleeding?

VETIGEL® stops moderate to severe bleeding in less than five seconds. In hard-to-reach places, light digital pressure can be used with a wet gauze to ensure contact with the bleeding site. Digital pressure may be used when dealing with severe bleeding cases. For more information, please consult IFU.

## Why VETIGEL® versus other hemostasis products?

When seconds count, your time is too valuable waiting for bleeding to stop. VETIGEL® lets you stop bleeding confidently, allowing you to focus on your patients and clinic's needs faster. Scan the QR code to learn more about VETIGEL®.



Learn more at  
[vetigel.com](http://vetigel.com)

Ref: PM-V0056-R01



## Storage Requirements

VETIGEL® should be stored dry and maintained between 2-25 C upon receipt. VETIGEL® should not be stored in direct sunlight.

## Is VETIGEL® Sterile?

Yes. VETIGEL® is terminally sterilized via gamma irradiation.

**SINCE 1981**  
**PENN**  
**Veterinary Supply, Inc.**  
Phone: 800.233.0210  
[www.pennvet.com](http://www.pennvet.com)

1-800-VETIGEL

[vetigel.com](http://vetigel.com)