




Innovation in Veterinary Styptics Stop the Bleeding Quickly and Easily


Nail trimming is an everyday event and challenge for all veterinary clinics. Accidentally trimming into the nail 'quick' with subsequent hemorrhaging is a daily occurrence. Pigmented nails or long nails have a variable length of quick making bleeding even more likely. It is important to stop this hemorrhage quickly and painlessly for the benefit of the patient and its owner.

I was lucky enough to try the HemaStop. It is the coolest little invention ever! It comes in a little compact tube that looks like Chapstick. Its easy to use and no mess.

Heidi Ward; DVM, ACVIM



Nail trim bleeding



Blot with gauze



Apply HemaStop

The **HEMASTOP™** stick is comprised of an astringent/styptic compound that is incorporated into a wax petrolatum substrate. The wax serves to seal the nail as well as hold the styptic in place. The wax substrate maintains **HEMASTOP™** in a tidy, usable form with none of the problems associated with the powdered styptics: no more rock-hard clumps or stained mess during application.



Packaged in Chapstick style and size screw tube.

In our practice we see many exotics, wild life, and of course dogs and cats. We have found the HemaStop product to be the best adjunct to nail trimming in all species. This efficacious product is clean with no messy powder and is so easy to use.

Jack Landess; DVM, Nokomis, FL



Phone - 800.233.0210

www.pennvet.com



Bird nail trim



Apply HemaStop



TRIMMING NAILS AND BEAKS

When confronted with a bleeding nail, one should hold the paw firmly, acting as a tourniquet, while blotting away blood. Then brush the **HEMASTOP™** into the cut surface. **One application is usually sufficient** but with very short trims **the process may have to be repeated** to stop the bleeding. For best results, use **HEMASTOP™** like putty, pressing it into the bleeding wound. Residual blood on the surface of the product will not affect its future usability, however it is easy to clean off by wiping it with a paper towel. This same application technique works when trimming beaks.

Some Interesting Facts About HEMASTOP™

- All the ingredients in **HEMASTOP™** are active ingredients. There are no artificial colors or fillers.
- Color variations in **HEMASTOP™** are normal and arise from small temperature fluctuations during the manufacturing process. This does not affect the efficacy of **HEMASTOP™**
- The inventor of the product is a practicing veterinarian with a half-century's experience in the profession. In his youth he was a boxer where he observed the "cut men" in the ring treating their fighters by using some of the principles that he later incorporated in the development of **HEMASTOP™**

When we have a toenail that has been clipped too close, we reach for HemaStop. It is efficient and there is no powdery mess to clean up afterward.

Chetika Maus; DVM, CVA
Sarasota, FL

Product # **J1257**
3pk
Manufactured for
Jorgensen Laboratories
Loveland, CO 80538
www.JorVet.com

Made in U.S.A.
Patent Pending

