

Stop chronic ear infections before they start.



Otter eQ: Quorum Sensing Technology

The Otter® from QSM Diagnostics is an advanced instrument using Quorum Sensing Technology designed to help veterinarians diagnose and monitor bacterial infections in dogs and cats.

Point-of-Care Testing


- **Diagnose the infection correctly the first time**
- **Provides fast results** - Results in as little as 5 minutes.
- **Helps provide better antimicrobial stewardship**
- **Additional revenue for your practice** - The Otter is a cost-effective solution for culture ID and can be bundled with your cytology for additional revenue.

Future cartridges:

- Multiplexed (7 common bacterial species)
- Fungal (3 common yeast species)
- Parasite (2 protozoa and 3 worms)




Otter eQ in action



QSM Diagnostics conducted a study to examine the accuracy of identifying Pseudomonas on cytology.

- 304 ear swabs collected from 9 different hospitals and specialty clinics
- First analyzed by cytology by LVTs or Veterinarians
- Sent to IDEXX BioAnalytics Lab for Bacterial Culture ID



SINCE 1981 PENN
Veterinary Supply, Inc.
 Phone: 800.233.0210
www.pennvet.com


Cytology Exam Result vs. Bacterial Culture ID Result

	YES Rod-shaped bacteria identified on culture	NO Rod-shaped bacteria identified on culture
Rods SEEN on cytology	199	59*
Rods NOT SEEN on cytology	14	32

*Includes 9 samples where no bacteria grew in culture

Accuracy, sensitivity and specificity of humans identifying rod-shaped bacteria during cytology exam

Accuracy	76%
Sensitivity	93%
Specificity	35%



Species of Rods on Ear Swab Culture after Cytology (N=258)

- Pseudomonas Only
- Pseudomonas with Other Rod Bacteria
- Other Rod Bacteria Only

How many of the rods “seen” on cytology are actually Pseudomonas?

- Cytology results marked rods on 258 of the 304 samples
- Only **29%** of these samples were positive for Pseudomonas in culture

