

Easy and accurate measurement of endotracheal tube cuff pressure.


Tru-Cuff

The affordable inflation syringe with a cuff pressure indicator

Endotracheal tube cuff pressure has typically never been monitored or measured in veterinary anesthesia. Most veterinarians and veterinary technicians are not aware of the suggested correct cuff pressure. Numerous studies in human hospitals have shown that over 60% of cuffs are over inflated with an average of 2-3 times the recommended guideline. Most clinicians go by a sense of feel or look at the one-way valve air pouch. Now they can choose the right product for this long overlooked, but important aspect of anesthesia.

Pressure Range:

Red Zone 40-60 cmH²O
Green Zone 18-26 +/-2 cmH²O




Why measure? The guideline for correct cuff pressure is 20-30cm H₂O. Pressure above 30cm H₂O will exceed the surrounding tracheal blood supply pressure and constrict proper blood flow. Prolonged exposure to high pressure will cause post operative coughing and soreness.

At pressure below 20cm H₂O, the under inflated cuff will allow aspiration of oral fluids especially during a dental procedure.

Instructions for use:

1. Attach the Tru-Cuff syringe to the endotracheal tube inflation valve.
2. Inflate until the black line is in the green pressure



J1404 3pk

LIT1404