

Matrix™ Central Scavenger

Waste Anesthetic Gas Disposal System



The Matrix™ Central Scavenger Scavenging Simplified.

The Matrix™ Central Scavenger...one less thing to worry about in your clinic.



Keep your mind free to focus on caring for your patients, not on scavenging excess anesthetic gas. The Matrix™ Central Scavenger safely removes waste gas from your clinic without the cost or hassle of charcoal canister replacement.

- Replace charcoal canisters forever with an active gas scavenging system that quietly and consistently keeps you safe
- Scavenge waste gas from up to 4 anesthesia machines
- Easily expand the system to cover up to 12 anesthesia machines without significant expense
- Keep your clinic neat and tidy without bulky interface devices attached to the anesthesia machines
- Simple installation process and flexible design ensures your scavenging system will be successfully installed
- Unparalleled service and support to give you peace of mind

Supported by Midmark

The Matrix™ Central Scavenger gives you unsurpassed reliability, backed by our trained service technicians and independent service network. You will have peace of mind knowing that over 90 years of manufacturing experience goes into each product we build. We have a veterinary design and engineering team that is rich in experience and focused on a customer-driven philosophy.

Specifications

Overall depth: 6" (15.24 cm)

Overall width: 13.5" (34.29 cm)

Overall height: 14" (35.56 cm)

Shipping weight: 28 lbs (12.7 kg)

Negative Pressure (vacuum) at zero flow – 0 to 0.5 cm H₂O

Occluded System (normal operation) - < 10 cm H₂O

Normal Operation, oxygen flush at 75L/min for 5 seconds - < 10 cm H₂O

Electrical Ratings

Voltage: 115 VAC +/- 10%

Model Numbers

40162800 Matrix™ Central Scavenger

Accessories

20236000 Expansion Kit (for installations with more than 4 machines)

20222500 Auto-Close Fitting Accessory, (19 mm)



Become of fan of Midmark Animal Health

Midmark is an ISO 13485 and ISO 9001 Certified Company. • 1-800-MIDMARK • 1-937-526-3662 • midmark.com
© 2009 Midmark Corporation • Products subject to improvement changes without notice • Litho in U.S.A. • 007-0819-00 Rev. A1 (10/09)

