



~~Tying
Gauze
Roll~~

EZ Tie

Endotracheal Tube Securing Device

So simple, you'll wonder why you did not think of it!!



Maintaining an endotracheal tube in the correct position is a critical component of anesthesia or surgery. Any shift or displacement of an endotube during anesthesia can be life threatening. A common practice is to use a length of cotton gauze roll to first tie around the tube then around the animal's muzzle. Roll gauze can come untied or be difficult to quickly untie.

*"Wow! I have been a technician for 20 years and intubated thousands of animals and now I'll never have to cut another piece of gauze or tie another knot. **EZ Tie** is ready to use every time and it's reusable!"* Tabitha Landow, Veterinary Technician, Ossining, NY

The **EZ Tie** was designed to be quick and easy to hold any endotracheal or anesthesia mask safely in place. The thin but strong plastic tubing has 2 sliding clasps. These clasps readily move along the tubing but hold steady under tension.

The first clasp is used to tighten around the ET tube or mask. The second slides along to tighten around the patient's muzzle. It also works extremely well with challenging brachycephalic breeds such as pugs and felines.






Patent Pending

J1202 EZ Tie pack of three

Will work with most anesthesia masks!



"I used the EZ Tie today and was impressed how easy it was to use for both cats and dogs and how secure it was with use on endotracheal tubes and anesthesia masks." Yvonne Szacki, DVM, Brooklyn, NY



The EZ Tie is:

- Simple to apply
- Reusable and durable
- Light weight
- Can hang in induction area ready for use
- Inexpensive
- One size fits all
- Watch short video on www.JorVet.com
- Not autoclavable

Jorgensen Laboratories, Inc.

Loveland, CO 80538

(800) 525-5614 Info@JorVet.com www.JorVet.com