

Eaze-Off™

Eaze-Off™ temporarily denatures the tack / adhesive properties of the glues used on high tack fixation, support, and compression bandages.

Allows easy and painless removal of the bandage or dressing minimizing trauma to the healing process and the animal. The tack of the glue is only temporarily suspended as Eaze-Off™ completely evaporates and leaves no residue or film barrier preventing the re-adhesion of the next dressing or new appliance to the animal.

CODE	DESCRIPTION	PACK
FE00050A	EAZE-OFF™ 1.69fl oz (50ml)	1
FE00200A	EAZE-OFF™ 6.76fl oz (200ml)	1

Features & Benefits



Dries in seconds allowing a replacement dressing to be applied immediately!



Scan to view an Eaze Off demonstration video.



Also works upside down!



Uses CFC Free Propellant.



STEP 1



STEP 2

Eaze-Off™ is totally harmless to the animal as Eaze-Off™ completely evaporates. Once evaporated from the bandage/dressing Eaze-Off™ becomes inactive and must be re-applied.



STEP 3



STEP 4



MILLPLEDGE VETERINARY

Tel: +1 (310) 215-0400 • Fax: +1 (310) 215-0401

5310 Derry Ave, Ste S&T • Agoura Hills • CA, 91301 USA

www.millpledge.com