

# 14687-28



Primary Plum Set<sup>™</sup>, Clave<sup>™</sup> Secondary Port, Clave Y-Site, Secure Lock, 103 Inch

## **Technical Specifications**

- >>Length: 103 in
- > Drop Former: 15 drops/mL
- > Approximate Priming Volume: 19 mL
- > Distal Length: 75 in
- > Approximate Distal Priming Volume: 13 mL
- > Tubing Material: Polyvinyl Chloride (PVC), non-DEHP

### Cassette Advantages

- > Air Trap can capture up to 1 mL of air, helping to minimize nuisance alarms
- Secondary port allows to easily remove air from the cassette without having to disconnect the IV Line from the patient, using the backpriming feature
- > Infusion of two compatible medications with a single dedicated set
- > Secondary port can also be used for syringe delivery

## Safety

- > Clave Needle Free connector in all ports of access help to minimize the infection risk for the patient
- > Flow regulator closes automatically if the pump door is opened to help prevent free-flow
- > Priming stop cap helps to prevent medication spills during priming
- > Components are not made with DEHP or natural rubber latex

#### Packaging

- > Case Quantity: 48 units per case
- > Shelf-life: 3 years from the date of production

#### **Replacement Requirements**

- > Plum Sets can be used up to 96 hours
- > Replace according to hospital policies

