



# ELDRIDGE ANESTHESIA MASK

A NEW INNOVATIVE MASK FOR OXYGEN THERAPY AND INDUCED GAS ANESTHESIA

The Eldridge anesthesia mask addresses a common problem. A standard round mask wants to roll away from a patient on a flat table. Towels are often used to brace the mask into position. The Eldridge mask is flat on one side to prevent this “rolling” and safely keep the mask in position.

## COMMON USES

- ▶ Induce gas anesthesia and sometimes maintain anesthesia in patients
- ▶ Oxygen therapy for severe respiratory distress

## FEATURES

- ▶ Clear, virtually unbreakable plastic cone
  - ▶ Allows visualization of patient’s muzzle area
- ▶ Stretchable black rubber diaphragm for snug fit
- ▶ 15mm OD end readily connects to any anesthesia circuit



J0355 Series

SKU	Animal	Mask Opening	Replacement Diaphragms SKU
J0355	Feline	1” diameter	J0355D #1
J0355A	Canine	1.75” diameter	J0355AD #2
J0355B	XL Canine	3” diameter	J0355BD #3
J0355S	Set of All 3		



Phone - 800.233.0210

www.pennvet.com

Jorgensen Laboratories, Inc.

800.525.5614 | fax 970.663.5042

www.jorvet.com | info@jorvet.com