





Revolutionary hydrophobic litter that keeps the cat's urine on top!



Noninvasive



100% Natural



Easy to Use



Cats Love it



Kit4Cat is the name of our patented hydrophobic litter that keeps cat urine on top, making sample collection stress-free! The "litter" smells and feels like regular sand with a natural non-toxic repellent coating. It allows veterinarians, technicians and cat parents to easily passively collect a urine sample without resorting to painful invasive procedures like Cystocentesis. The kit comes with everything you need to collect a urine sample.

Kit4Cat is approved by veterinarians around the world, having gone through extensive validation testing. Since 2010, Kit4Cat has been used in clinics and hospitals in 44 countries worldwide. It is sold by every major veterinary distributor. Because this method is easy to facilitate and the acceptance rate is so high, Kit4Cat is regarded as the most humane and ideal method for noninvasive urine collection.



Manufactured by: Coastline Global, Inc. 827 Lincoln Ave, Suite B11, West Chester, PA 19380 Tel: 888-999-3050 | info@coastlineglobal.com





Benefits for Cat Owners:

- Noninvasive and easy to use at home, keeping cats stress-free
- No hospitalization required for urine sampling
- Sand is biodegradable, so disposal is easy
- Cats love it because it feels like regular sand
- Odorless & non-toxic
- Cat will remain stress-free and happy
- Sedation and hospitalization not required
- Sample can be collected at home Great for monitoring diabetic cats
- Perfect for cats with chronic kidney disease
- · Available from US and Canadian veterinarians

Benefits for Cats:

- 100% natural sand
- Painless and stress-free urine collection
- Non-invasive
- Non-toxic
- Ododess





Benefits for Veterinarians & Technicians:

- No more forcing cats to urinate with a mess to clean up
- No more cat bites or scratches

- Recommended by top veterinarians
- Another quality product for you to sell
- Available from top US & Canadian distributor



Everybody loves Kit4Cat

"We have started to use the Kit4Cat sand in place of the Nosorb. Clients and their cats like it more. We love the product." Dr. Cynthia Wojcicki, Merrimac Valley Animal Hospital.

"We thought the kit was great."

Karen Campbell, Veterinary Technician, Smithtown Valley Animal Hospital

"4 days with Nosorb litter, cat kept urinating on paper or blanket. When sand came in, put in cage and overnight cat used box, able to obtain free catch sample."

Andrea Stoops, Vet Nurse, Ayer Animal Medical Center

In the box

3x 11oz bags of Kit4Cat hydrophobic sand 3x disposable pipettes for urine sample collection 3x disposable vials with caps for holding urine sample



Pour the litter Pour all hydrophobic litter in a clean and empty litter box.



Fill the vial Fill the vial at least 3/4 way using the pipette.



Wait for your cat Allow your cat to urinate freely on the hydrophobic litter.



Store the sample Close the vial with the cap and keep in refrigerator for up to 48h.



Collect sample Easily collect urine from the top of the litter with the pipette.



Further analysis





Vet Approved