

PinnaMend™

The Original Sutureless Aural Hematoma Repair Kit



Repair done in less than 10 minutes with
no bandaging required.



"Partnering with you to bring New Technology, Quality, Reliability and Service
Above and Beyond"


- Repair done in less than 10 minutes
- No bandaging required
- The patients ear is secured to the top of the patients head with the PinnaButton
- The button method reduces swelling leading to faster healing times and keeps ear canal open
- Fenestrations on pad allow for skin to breath System can be removed 2 weeks after procedure
- Removal of PinnaMend is quick & painless
- Low rate of re-occurrence
- Low rate of cauliflower ear as a result of uniform scar tissue formation



PinnaMend Starter Kit

Everything you need for your first procedure:

- 2 Silicone Ear Pads
- Attachment Buttons
- 20 Quick Lock Clips
- 40 Quick Lock Rings
- A 14 Gauge Needle (Stylet)
- Surgical Procedure Guide



PinnaMend Refill Kit

- Sterile packs of the clips rings
- Needle to do 4 additional procedures
- Pads & Button are autoclaved and reused

