

SIZE CHART

ITEM	DESCRIPTION	SIZE	MAX CIRC.
PVSFF8	PennChoice FastFit 8cm	X-Small	10.5"
PVSFF12	PennChoice FastFit 12cm	Small	12.5"
PVSFF17	PennChoice FastFit 17cm	Medium	16"
PVSFF23	PennChoice FastFlt 23cm	Large	19.5"
PVSFF30	PennChoice FastFit 30cm	X-Large	24.5"

FITTING: Prior to applying Velcro®, place collar on dog's neck to measure for Velcro® strip placement. Then, place Velcro® on opposing sides of the collar. Once applied, Velcro® is not adjustable.

