



KVP Rehabilitation

# Balto® Hock

## Hock Brace

Addressing hock pain, lameness and neuromuscular problems, the Balto® Hock also aides:

- Achilles injuries and orthopedic conditions
- Neurological conditions
- Pre/Post-surgical immobilization
- Contralateral support
- Protection and supporting bandages

The brace can even be used as an alternative to surgery if the patient is not an ideal candidate and can also be used without splints for conditions such as arthritis and arthrosis. Compression helps increase blood pressure in the area, raising the temperature of the immediate area and helping to reduce pain.

ITEM #	SIZE	WEIGHT	A	B	C
BTHOCKXXS	XXS	5 - 10 lbs	3.75"	3"	1.75" - 2.75"
BTHOCKXS	XS	10 - 20 lbs	4.25"	3.25"	2.25" - 3"
BTHOCKS	S	20 - 55 lbs	5"	4"	3" - 4.25"
BTHOCKM	M	55 - 75 lbs	5.75"	4.75"	4.25" - 5"
BTHOCKL	L	75 - 110 lbs	6"	5"	5" - 5.75"
BTHOCKXL	XL	over 110 lbs	6.25"	5.5"	5.75" - 6.25"

**MEASUREMENT A:** Circumference above the hock

**MEASUREMENT B:** Circumference below the hock

**MEASUREMENT C:** Height of the brace

