

Which Medipaw® Boot Should I Use?

Basic Boot, X Boot and Slim Boot Comparison Chart



Medipaw® Basic Boot

Used for: Covering soft bandages

Duration: 2 weeks or less

Goal: To provide comfort and protection from water, snow, dirt, mud and rain.

Boot Details:

- Color-coded sizing – available in XS, XS2, SH, S, M, L, XL, RJ
- Adjustable elastic Velcro® straps
- Light, flexible, & durable rubber sole
- Upper drawstring closure
- Waterproof/breathable nylon shell
- Machine washable/hang dry only
- Seam-sealed
- Scuff guard



Medipaw® X Boot

Used for: Covering thicker bandages, casts or splints

Duration: More than 2 weeks

Goal: To accommodate thicker dressings or patients with difficult gaits on difficult terrain, for longer periods of time.

Boot Details:

- Color-coded sizing – available in XS, XS2, SH, S, M, L, XL, RJ
- Scuff guard
- Adjustable elastic Velcro® straps
- Light, flexible, & durable rubber sole
- Upper drawstring closure
- Waterproof/breathable nylon shell
- Machine washable/hang dry only
- Seam-sealed
- Scuff guard



Medipaw® Slim Boot

Used for: Protecting wounds without bandages

Duration: 2 weeks or less

Material: Antimicrobial bamboo and silver

Goal: To reduce bacterial threats and improve healing. Helps with: burns, degloving, excessive licking and irritation, hot spots, lick granulomas, punctures.

Boot Details:

- Color-coded sizing – available in XS, XS2, SH, S, M, L, XL
- Machine washable/hang dry only
- Durable, non-skid, rubber sole
- Upper drawstring closure
- Adjustable elastic Velcro® straps
- Reflective strap
- Breathable nylon shell
- Scuff guard

We've got you covered!

US: (800) 245-3413 | The Netherlands: +31 646 3344 45
info@medivetproducts.com | www.medivetproducts.com



Which Medipaw® Boot Should I Use?

Basic Boot, X Boot and Slim Boot Comparison Chart



	Medipaw® Basic Boot	Medipaw® X Boot	Medipaw® Slim Boot
Sizes Available:	XS, XS2, SH, S, M, L, XL, RJ	XS, XS2, SH, S, M, L, XL, RJ	XS, S, M, L, XL
Used to cover:	Soft bandages	Thicker bandages, casts or splints	Open wounds
Can be used for:	2 weeks or less	More than 2 weeks	2 weeks or less
Adjustable Velcro® Straps	✓	✓	✓
Upper Drawstring Closure	✓	✓	✓
Machine Washable	✓	✓	✓
Seam-sealed	✓	✓	✗
Reflective Strap	✗	✗	✓
Scuff Guard	✓	✓	✓
Waterproof/breathable	✓	✓	✗

We've got you covered!

US: (800) 245-3413 | The Netherlands: +31 646 3344 45
info@medivetproducts.com | www.medivetproducts.com