

HEMOSTEM™

HEMOSTATIC

Dental Cubes with carboxycellulose



FOR USE IN DENTAL PROCEDURES

Control protracted postextractional capillary bleeding.

EASY INSERTION



CAN BE
SUTURED IN
PLACE

Absorbed by the body within 8 days- no need for it to be removed. Additional cubes can be used at the same site to control bleeding quicker or for covering larger areas or deeper and wider sockets.

CONTAINS:




ANTIMICROBIAL EFFECT



PRESENTATION

5 CUBES



STERILE

CODE	QUANTITY	PACK
HP00510A	5 CUBES	1

**STOP
bleeding**

HEMOSTEM™ is used for the control of bleeding, as capillary haemorrhaging, and serum exudating from damaged tissue wherever conventional methods of local haemostasis is difficult or impractical



HemoStem™ is not suitable for the control of arterial or venous pressure bleeding where appropriate prior 'tying off' or cautery of blood vessels must be carried out. It may be used to control small arteriolar or venous bleeding in combination with application of short time digital pressure or pressure packing.

HEMOSTATIC

Gauze with carboxycellulose



FOR USE DURING SURGERY, OPEN CAVITY, LARGER BLEEDING AREAS AND SURGICAL PROCEDURES



Delicately woven hemostatic gauze for wrapping around anastomoses such as those of arterial and venous vessels and the digestive tract.

STOPS BLEEDING FAST



Use sparingly. One layer is normally required for wound surface application. For packing into soft or hard tissue wound cavities, use only sufficient to adequately cover inner cavity surfaces.



PRESENTATION

2" x 2¾"

STERILE

CODE	SIZE	PACK
HP00505A	2" x 2¾"	1

DIRECTIONS: Cut open peel pouch and carefully remove gauze with forceps. Cut to size as required.



**MILLPLEDGE
VETERINARY**

Tel: +1 (310) 215 0400 • Web: • Email: nasales@millpledge.com

www.millpledge.com



53 Industrial Circle • Lancaster • PA 17601 • 1-800-233-0210 • www.pennvet.com