

SIZE CHART

ITEM	DESCRIPTION	CIRCUMFERENCE	DEPTH
PVSBTC175	PennChoice Clear 7.5cm	9" 10.25"	3"
PVSBTC210	PennChoice Clear 10cm	9.25"- 10.5"	4"
PVSBTC312	PennChoice Clear 12.5cm	11"– 12.25"	5"
PVSBTC415	PennChoice Clear 15cm	13"– 16"	6"
PVSBTC520	PennChoice Clear 20cm	13.25"- 16.25"	8"
PVSBTC625	PennChoice Clear 25cm	15.5"– 19"	10"
PVSBTC730	PennChoice Clear 30cm	18"- 22.25"	12"
PVSBTC840	PennChoice Clear 40cm	25"- 21"	16"

NECK CIRCUMFERENCE: Wrap a soft measuring tape around the pet's neck, where their collar would sit.

COLLAR DEPTH: Measure the distance from the pet's nose to the base of their neck, just before the shoulder blade.

