


Get your patients on the mend faster than sutures.

Closing wounds with sutures can be time-consuming and costly. 3M™ Vetbond™ Tissue Adhesive lets you seal wounds quickly, accurately and conveniently. So your patients can get back to doing what they love faster.



3M™ Vetbond™ Tissue Adhesive

3M™ Vetbond™ Tissue Adhesive quickly closes minor wounds, often eliminating the need for sutures or bandages. And now it's paired with the next generation of adhesive packaging, making it faster, safer and more convenient to use than ever.

- Economical 3 mL bottle includes four “reusable” applicator tips
- Blue dye for easy-to-see drop application
- Unique nozzle seals itself with a proprietary pin cap and eliminates the need for a cutting blade to open
- Redesigned bottle and cap keeps the nozzle clear and allows unlimited opening and resealing, maximizing adhesive use
- Yellow pipettes ensure easy visibility when placed on a surgical table or white surface
- Ideal for a range of procedures, including cat declaws, dew claw removal, closing minor lacerations, tooth extractions and oral surgery



Ordering information

Product Number	Product Name	Size	Packaging
1469SB	3M™ Vetbond™ Tissue Adhesive (For Veterinary Use Only)	3 mL	12 bottles/case

Learn more about 3M's innovative solutions for animal care at 3M.com/animalcare



Animal Care

3M Center, Building 223-5N-03
St. Paul, MN 55144 USA

Phone 1-800-848-0829
Web 3M.com/animalcare
3M.com

3M and Vetbond are trademarks of 3M.
Please recycle. Printed in USA.
© 3M 2016. All rights reserved.
70-2011-5823-8



53 Industrial Circle • Lancaster • PA 17601 • 1-800-233-0210 • www.pennvet.com