


Clinical signs of canine hypothyroidism

Canine hypothyroidism


Lethargy
Mental dullness
Weigh gain
Inactivity
Weakness
Anemia
Endocrine alopecia
Behavioral abnormalities

Clinical signs are variable and depend in part on the age of the dog at the time a deficiency in thyroid hormone develops.




Test procedure


1 SAMPLE
100µl



2 MIX
10 times



3 LOAD
100µl



Features



Sample
Serum 100 µl



Measurement Range
0.25-5.00 mg/L



Test time
15 minutes

Ordering Information

Product No.	Product Name	Product Type	Packing Unit
VCF118DC	Vcheck cTSH	Device	5 Tests/Kit



To learn more and place your order

Minneapolis Minnesota

✉ customerservice@bionote.com ☎ 800-727-5169