

# ULTRoS



| Unit                        | UltrOs X5  | UltrOs MCX5  | UltrOs LX5   | UltrOs CX5   |
|-----------------------------|--|--|--|--|
| Configuration               | Convex 50mm and Linear 25mm, Wireless                                      | MicroConvex, Wireless  | Linear 40mm Width  | Convex 60mm  |
| Frequencies                 | Convex probe 3.5/5.0 MHz, Linear probe 7.5/10 MHz                          | 6.5/8.0 MHz  | 7.5/10 MHz   | 3.5/5.0 MHz  |
| Size and/or useful scanning | 5-150 Lbs. 0.5 - 25 cm's   | 10-80 lbs. 1-10 cm's   | 0.5- 6 cm's  | 30-220 lbs. 9-30 cm's  |
| Scan Modes                  | B, B/M, Color, PW  | B, B/M, Color, PW  | B, B/M, Color, PW  | B, B/M, Color, PW  |
| Size                        | 4.75 x 2.0 x .75 in.   | 5.0 x 2.2 x 0.8 in.  | 6.4 x 2.2 x 0.8 in.  | 4.5 x 2.8 x 0.8 in.  |
| Weight                      | 4 oz.  | 3.2 oz.  | 3.5 oz.  | 4.4 oz.  |
| Image Adjust                | Gain, Depth, DR, Focus, SRI, Harmonic                                      | Gain, Depth, DR, Focus, SRI, Harmonic                            | Gain, Depth, DR, Focus, SRI, Harmonic                            | Gain, Depth, DR, Focus, SRI, Harmonic                            |
| Measure                     | Length, Area, Velocity, HR, S/D, Obstetrics, Auto Blood Vessel Measurement |  |  |  |
| Puncture Assist Function    | guide line, out-of-plane puncture guide line, needle enhance               |  |  |  |
| Battery                     | Built-in lithium battery, wireless or USB charge 1.5h, work 1.5h           | Built-in lithium battery, wireless or USB charge 1.5h, work 1.5h | Built-in lithium battery, wireless or USB charge 1.5h, work 1.5h | Built-in lithium battery, wireless or USB charge 1.5h, work 1.5h |

Image

