



**COOLRENEWAL**  
www.cool-renewal.com

*Providing a Low Cost,  
High Value  
Cryosurgical Solution  
to Veterinarians  
around the World*



### Benefits

- Minimal Discomfort
- Self-Anesthetizing
- 2 Minute Procedure
- No Suturing
- No Evaporation
- No Expiration
- Portable
- Low Cost Treatment
- 3 Disposable Applicator Options
- Multiple Applicator Sizes for Various Lesion Sizes and Locations

### Commonly Treated Skin Lesions

**Papillomas**



**Sebaceous Adenoma**



**Skin Tag**




**Distichia**



**Histiocytoma**



**Sarcoid**



Other treatable lesions include, but are not limited to:

**Follicular Cysts, Granulomas, Fibroids, & more!**

Cool Renewal™ is a convenient and cost effective alternative to traditional liquid nitrogen methods and equipment, as well as other expensive portable cryosurgery systems. Cool Renewal™ was developed to fill the gaps of its competitor products, offering veterinarians a wide variety of applicators in several sizes, new applicators unique to our cryosurgical system, and the lowest cost per treatment available from a portable cryosurgical system.

## Cool Renewal™ Complete Cryosurgical Kit

Item # J1401



- ❖ (2) 170mL Canisters of Cryogen  
~ 130 Freezes\*
  - ❖ (2) Extender Tubes
  - ❖ Assorted Pack of 60 Disposable, Double Ended Foam Tipped Applicators  
(20 ea, Small, Medium, Large/XL)
  - ❖ Assorted Pack of 50 Disposable Plastic Isolation Funnels  
(10 ea- 5mm, 8mm, 12mm, 14mm, 16mm)
  - ❖ (20) Disposable Skin Tag Tweezers
  - ❖ Instructions for Use
  - ❖ Patient Instructions Tear Pad
  - ❖ Applicator Practice Pad
  - ❖ Magnetic Packaging for Component Storage
- \*Average Number of Freezes based on treatment of 5mm lesions

## Cool Renewal™ Replacement Canister

Item # J1401B



- ❖ Single Replacement Canister-  
170ml (65 Freezes)
- ❖ Extender Tube
- ❖ Instructions for Use
- ❖ Patient Instructions Tear Pad

### Foam Tipped Applicators



- ❖ Small-Item# J1401D1
- ❖ Medium- Item# J1401D2
- ❖ Large/ XL-Item# J1401D3
- ❖ Assorted- Item# J1401D4  
\*60 per pack

### Skin Tag Tweezers



- ❖ Item # J1401A  
(qty 20)

### Isolation Funnels



- ❖ 5mm- Item#J1401D5  
(qty 50)
- ❖ 8mm- Item#J1401D6  
(qty 50)
- ❖ 12mm- Item#J1401D7  
(qty 50)
- ❖ 14mm- Item#J1401D8  
(qty 50)
- ❖ 16mm- Item#J1401D9  
(qty 50)
- ❖ Assorted-  
**Item#J1401D10**  
(10 of each size)

All applicators are disposable. DO NOT ATTEMPT TO CLEAN OR USE ON MULTIPLE PATIENTS

Jorgensen Laboratories  
[www.JorVet.com](http://www.JorVet.com) info@JorVet.com

Please  
Contact  
Us At:

**PENN**  
Veterinary Supply, Inc.

53 Industrial Circle  
Lancaster, PA 17601

Phone: 1-800-233-0210  
Online: [www.pennvet.com](http://www.pennvet.com)