




DOPPLER MONITOR PROBE GEN II

HYPERTENSION

NO DISEASE CAN BE DIAGNOSED
IF YOU DON'T LOOK FOR IT!



FEATURES

- ▶ Lightweight, Hand-held Unit
- ▶ High Quality, Trigger-style Sphygmomanometer
- ▶ Rechargeable Lithium Ion Battery – Approx. 6 Hours of Use w/Full Charge
- ▶ Includes 5 Sizes of Blood Pressure Cuffs, Headphones, and Tube of Gel
- ▶ 4-minute Auto Shut Off
- ▶ Complete w/Hard-sided Carrying Case
- ▶ 1-year Limited Warranty on Doppler and Probe Only

SKU

- J0563N** Complete Unit w/8 mHz Flat Probe
- J0563ND1** 8 mHz Flat Probe
- J0563ND3A** Sphygmomanometer Only
- J0563ND10** Pencil Probe for Exotics
- J0563ND11** Esophageal Probe - Optional

PROBLEM

Hypertension (High Blood Pressure) is often overlooked during diagnosis, treatment, and routine physical exams in small animals. In addition, many hypotensive (Low Blood Pressure) patients, such as those suffering from shock, benefit greatly from Doppler measurement of blood pressure when knowing the pressure can mean effective and timely treatment.



SOLUTION

Combining Doppler with an appropriate size “cuff” for small animals, thus providing the veterinarian with the same, noninvasive method for testing blood pressure as is available in human medicine.



TREATMENT

Newer drugs are being used for treatment of hypertension. There is a need to monitor the patient’s response to treatment with a Doppler for precise readings.

DOPPLER MONITOR GEN II PROBE OPTIONS



SKU | J0563ND1
8 mHz Flat Probe

- ▶ Contoured for positioning on paw or tail
- ▶ Large beam with high sensitivity
- ▶ Robust compact design



SKU | J0563ND3A
Sphygmomanometer

- ▶ High-quality human medical grade
- ▶ Trigger-style release



SKU | J0563ND10
Pencil Probe

- ▶ Used to locate and monitor pulses in exotics and neonatal patients
- ▶ Aid in diagnosis of arterial disease



SKU | J0563ND11
Esophageal Probe

- ▶ Cardiac monitoring solution for intubated patients
- ▶ A must have when monitoring critical anesthetized patients