

Medium support, highly conforming knitted, compression and support bandage. Specifically designed for multi-purpose animal bandaging applications.

KNIT-FIX™ A light-weight open-knitted cellular structure for good ventilation. Providing excellent conforming, non slip and non fray properties for moderate support and holding dressings in position.

Greater length and width ways stretch gives more conformability and distribution of the bandage. This offers better sub-bandage pressure and limb coverage around awkward areas such as elbow and hock joints, giving additional comfort to your patient.

2ND LAYER

Perfect for use as the second layer in the three layer bandaging technique.

EXTRA STRETCH GIVES **EXTRA LENGTH**

Over 40% more usable stretched length when compared to other similar premium alternatives.

> A MEDIUM WEIGHT BANDAGE Ideal for

More flexible and more versatile for the variables of multi species practices.

MEDIUM multi purpose use.

FLEXIBLE & TRULY CONFORMS

More conformability and even distribution of bandage offering superior support and limb coverage.

No slipping, no sagging.

FEATURES & BENEFITS

NO LOSS OF SUPPORT FOR **EXTENDED** PERIODS OF USE

Superb for lightweight protection to compression and immobilization bandage techniques.

SUPERIOR STRUCTURE

Gives Warmth, Comfort and less risk of slipping, providing greater security.



46% more usable length than Conform. More bandage for your money!



WRAPPED INDIVIDUALLY

Bandages stay cleaner in storage.

Lightly elastic with some sideways extensibility and up to 70% lengthways extensibility.

Compared to any woven materials, the unique knitted structure provides a veterinary bandage that is very strong.

Available in a wide range of widths. Choose the correct size bandage to fit the animal

CODE	DESCRIPTION	PACK
DK01025A	KNIT-FIX™ 1" x 2.25yds	10
DK01050A	KNIT-FIX [™] 2" x 4.4yds	10
DK01075A	KNIT-FIX" 3" x 4.4yds	10
DK01100A	KNIT-FIX" 4" x 4.4yds	10
DK01150A	KNIT-FIX" 6" x 4.4yds	10

Distance of relaxed lengths and new distances of stretched lengths:

Knit-Fix™ relaxed length 108" / 3yds	~~~~~
Knit-Fix™ stretched length 158" / 4.4yds)
Knit-Fix™ recovery length after Stretch = 100% recovery	
Conform ^e relaxed length 75" / 2yds	
Conform® stretched length 108" / 3yds)
Conform® recovery length after Stretch = only 74% recovery	├



MILLPLEDGEVETERINARY

5310 Derry Ave, Ste S&T, Agoura Hills, CA, 91301 USA Tel: +1 (310) 215-0400 Fax: +1 (310) 215-0401 Email: nasales@millpledge.com

