



PennChoice Snap™

- Provides full body coverage
- Simple tabbed assembly for an easy fit
- Clear, lightweight, resilient material
- Curable plastic withstands cracking in cold weather
- Made in the USA

SIZE CHART

ITEM	DESCRIPTION	CIRCUMFERENCE	DEPTH
PVSKQ75	PennChoice Snap 7.5cm	9"– 10.25"	3"
PVSKQ10	PennChoice Snap 10cm	9.25"– 10.5"	4"
PVSKQ125	PennChoice Snap 12.5cm	11"– 12.25"	5"
PVSKQ15	PennChoice Snap 15cm	13"– 16"	6"
PVSKQ20	PennChoice Snap 20cm	13.25"– 16.25"	8"
PVSKQ25	PennChoice Snap 25cm	15.5"– 19"	10"
PVSKQ30	PennChoice Snap 30cm	18"– 22.25"	12"
PVSKQ35	PennChoice Snap 35cm	22"– 26"	14"
PVSKQ40	PennChoice Snap 40cm	25"– 31"	16"

NECK CIRCUMFERENCE: Wrap a soft measuring tape around the pet's neck, where their collar would sit.

COLLAR DEPTH: Measure the distance from the pet's nose to the base of their neck, just before the shoulder blade.