


14220-28

Primary Plum Y-Type Blood Set, 200 Micron Filter, Clave™ Secondary Port, Clave Y-Site, Non-Vented, Secure Lock, 110 Inch



Technical Specifications

- › Length: 110 in
- › Drop Former: 20 drops/mL
- › Approximate Priming Volume: 52 mL
- › Distal Length: 77 in
- › Approximate Distal Priming Volume: 16 mL
- › Tubing Material: Polyvinyl Chloride (PVC), non-DEHP

Cassette Advantages

- › Air Trap can capture up to 1 mL of air, helping to minimize nuisance alarms
- › Secondary port allows to easily remove air from the cassette without having to disconnect the IV Line from the patient, using the backpriming feature
- › Infusion of two compatible medications with a single dedicated set
- › Secondary port can also be used for syringe delivery

Safety

- › Clave Needle Free connector in all ports of access help to minimize the infection risk for the patient
- › Flow regulator closes automatically if the pump door is opened to help prevent free-flow
- › Priming stop cap helps to prevent medication spills during priming
- › Components are not made with DEHP or natural rubber latex

Component Highlights

- › T-Type configuration allows the connection of two different bags to the primary line
- › 200 micron filter in the drip chamber is designed to remove blood clots and harmful particles

Packaging

- › Case Quantity: 48 units per case
- › Shelf-life: 3 years from the date of production

Replacement Requirements

- › It is recommended to change sets within 4 hours after initiating infusion of blood products
- › Replace according to hospital policies