

# Vet Pump

## **Volumetric Infusion Pump**

With IV fluid warmer

J1060Q




- \*Use any standard IV set.
  Store up to 40 brands in
  memory
- \*IV fluid warmer
- \*Affordable-NEW, not refurbished
- \*Easy to read, bright LCD screen
- \*Easy to navigate menu options
- \*Lightweight and portable, < 5 lbs
- \*Programmable infusion

Veterinary Supply, Inc.

Phone: 800.233.0210 www.pennvet.com

### **Ease of Use of the JorVet Infusion pump:**



### **Components and functions:**





horizontal or vertical bars

- \*Runs on standard 110V AC or battery.
- \*Simple to navigate menu screen lets you set infusion parameters, as well as change speaker volume, occlusion pressures, air sensitivity and calibrate IV sets.
- \*IV fluid warmer- select warmer settings in menu and simply wrap the line around in the warming grooves.
- \*Stores past infusion records

### Specifications for the J1060Q:

- Easy to navigate menu screens
- Pumping mechanism: Peristaltic finger pump (NOT FOR BLOOD)
- Applicable IV set: 10, 15, 20 and 60 d/ml (with outer diameter 3.4 ~ 4.5mm)
  Can program up to 40 different brands of dedicated IV sets.
- IV fluid line warmer: 25-40°C (77-104° F) settings. IV tubing wraps around warmer through groves.
- Flow rate setting range: 1.0 -1200 ml hr (can set in 0.1ml increments up to 100 ml)
- Infusion rate options: Set at ml/hr, drops/ minute or volume/time
- Flow rate accuracy: +/- 5% (using calibrated standard IV set)
- Volume limit: 1-9999 mls
- KVO: Adjustable 1 ml/hr-20ml/hr
- Purge rate: 600 ml/hr
- Bolus mode
- **Programable**: Up to 8 different infusion segments
- Audio and visual alarms: Air-in-line, Door open, Upstream occlusion, Downstream occlusion, Infusion finished, No operation, Battery low, Pump error, Motor reverse, AC power lost, fluid warmer failure, flow rate error.
- Self-check during start up.
- **Power source**: AC100V ~ 240V, 50Hz/60Hz
- Internal battery: Rechargeable Ni-MH battery, DC9.6V, 2000mAh. Approx. 4 hrs at 25ml/hr
- Display: LCD with high brightness
- **Dimensions:** 13cm x 12cm x 19.5cm
- Weight: 2 kg
- IV pole mount: Adjustable for horizontal or vertical bars
- 1 year limited warranty



Jorgensen Laboratories, Inc. 1450 Van Buren Ave Loveland, Colorado 80538

Telephone: (800) 525-5614 <u>www.jorvet.com</u> Fax: (970) 663-5042 info@jorvet.com