



PennChoice Classic™

- Provides full body coverage
- Multiple adjustment points for a proper fit
- Secures to collar for additional support
- Curable plastic material withstands cracking in cold weather
- Made in the USA

SIZE CHART

ITEM	DESCRIPTION	CIRCUMFERENCE	DEPTH
PVSBT175	PennChoice Classic 7.5cm	9"– 10.25"	3"
PVSBT210	PennChoice Classic 10cm	9.25"– 10.5"	4"
PVSBT312	PennChoice Classic 12.5cm	11"– 12.25"	5"
PVSBT415	PennChoice Classic 15cm	13"– 16"	6"
PVSBT520	PennChoice Classic 20cm	13.25"– 16.25"	8"
PVSBT625	PennChoice Classic 25cm	15.5"– 19"	10"
PVSBT730	PennChoice Classic 30cm	18"– 22.25"	12"

NECK CIRCUMFERENCE: Wrap a soft measuring tape around the pet's neck, where their collar would sit.

COLLAR DEPTH: Measure the distance from the pet's nose to the base of their neck, just before the shoulder blade.