

800-640-7780 www.MicroVetDiagnostics.com





The InSight Urinalysis reader is a semiautomated urine test strip reader with touch screen and printer to be used in conjunction with InSight Urine Test Strips. Dip the strip into the urine and place it on the strip loader ....have your results in 60 seconds.

UPC results & Ratio with in 2 minutes.

Microalbumin and Creatinine results
& ratio within 2 minutes.

## **MS-11 Strips**

Urobilinogen, Blood, Glucose, pH, Bilirubin ,Protein, Ketone, Creatinine, Specific Gravity, Nitrate, Leucocytes

## **Protein/Creat Ratio**

## MS-2 Strips

Microalbumin & Creatinine
With Microalbumin/Creat Ratio

