


# Olympic Clean-Catch™

**A faster, easier way  
to collect clean urine specimens**



**CLEAN**  
disposable  
cups

### **Easy to Use – Saves Time**

The Olympic Clean-Catch makes urine collection faster and easier with dogs of all sizes – and you get clean specimens. With males, it is easy to position the Clean-Catch to collect clean mid-stream specimens. You can see what you are doing and make sure you get every drop. With females, the narrow, low profile cup slips easily under even small dogs.

### **Dog Friendly**

Clean-Catch saves valuable minutes with every urine collection. Trying to collect urine with a hand-held cup or ladle is slow and inefficient. Some dogs get nervous and won't urinate when you are standing too close to them or will stop urinating when you place a cup under them. But the Clean-Catch is dog friendly. It keeps you at a distance, so the dog is more likely to urinate, and it won't spook the dog when you put it under him or her.

### **No More Bending Down to Collect Samples**

When you use the Clean-Catch, you don't have to bend over or stoop down. Easy on your back. The long handle and flat cup give you visibility, so you get a good specimen every time. More sanitary and makes urine collecting less messy.

### **Clean-Packed Cups. No Cross Contamination.**

When collecting mid-stream urine samples it is extremely important to have a clean container to avoid possible cross contamination between patients (bacterial, cellular, artifactual). The plastic cups furnished with Clean-Catch are clean-packed and for one patient use only. This ensures more reliable lab results and improved patient care. The holder also keeps the cups from touching the ground. The cups have a handy pouring spout and are made of enviro-friendly plastic.

### **Eliminates Washing – Always Ready for Use**

With these inexpensive disposable cups, no time is wasted washing cups or ladles. Each disposable cup costs less than a few seconds of a tech's salary. And since no washing is required, the Clean-Catch collector is always ready for the next patient.


### **Pays for Itself**

The new Olympic Clean-Catch saves you time and money in three ways – you can collect specimens faster, it eliminates washing, and it is always ready for immediate re-use. By cutting costs, it will quickly pay for itself.

# Olympic Clean-Catch™

Clean-packed cups prevent cross-contamination.  
Single patient use assures more reliable lab results.  
Eliminates washing – saves time.

- **Easier – no bending – no stooping**
- **Dog friendly – saves time collecting urine**
- **Ideal for males and females**
- **Collect clean mid-stream samples**
- **Always ready for the next patient**
- **Saves time – saves money**



Clean, single use cups protect patients




Just push cup in

Lift cup out



Recyclable plastic  
Enviro-friendly

Narrow oval, less than  
1.4" high – ideal for small  
females and all males



Patent Pending



## Order From Your Distributor

Olympic Clean-Catch Handle Only ..... Model 50101

Olympic Clean-Catch Disposable Cups ..... Model 50105  
Box of 120