



Protexis[™] Latex Hydrogel

- Balances tactile sensitivity with protection, even when double-gloving
- Water-based hydrogel coating for easy donnability with wet or dry hands and enhances second-skin comfort of latex

Catalog no.	Size	Length	Thickness ¹			Material	Calan	Cuff	Protein	Qty/	Qty/
			Finger	Palm	Cuff	Material	Color	type	content ²	bx	cs
2D72LS55	5.5	11.5 in./ 292mm	9.8 mil/ 0.25mm	9.1 mil/ 0.23mm	7.1 mil/ 0.18mm		Translucent Yellow	Beaded/ Rolled	Less than 50 μg/dm²	50	200
2D72LS60	6					Natural rubber latex coated with acrylic hydrogel polymer coating					
2D72LS65	6.5										
2D72LS70	7	12 in./ 305mm									
2D72LS75	7.5										
2D72LS80	8										
2D72LS85	8.5										
2D72LS90	9										49/

See Appendix page 31 for complete testing standards



Phone - 800.233.0210

www.pennvet.com



SURGICAL GLOVES

^{2.} Protein content tested using ASTM D 5712